

PDO (Polydioxanone) Bi-Directional

ABSORBABLE

Precision Reverse Cutting

Length

	Color	cm	in	4-0	3-0	2-0	0	1	2
P-1  12 mm 3/8 circle	Violet	7 x 7	3 x 3		RA-1017Q				
PS-2  19 mm 3/8 circle	Violet	3 x 3	1 x 1		RA-1008Q	RA-1002Q			
		7 x 7	3 x 3	RA-1012Q-0	RA-1010Q	RA-1007Q			
		14 x 14	5 x 5		RA-1050Q-0	RA-1046Q			
		24 x 24	9 x 9		RA-1051Q-0	RA-1036Q-0			
PS-1  24 mm 3/8 circle	Violet	45 x 45	18 x 18				RA-1090Q		
PS  26 mm 3/8 circle	Violet	7 x 7	3 x 3			RA-1016Q	RA-1015Q		
		14 x 14	5 x 5		RA-1021Q	RA-1020Q-0	RA-1019Q-0		
		24 x 24	9 x 9		RA-1052Q-0	RA-1028Q-0	RA-1030Q-0		

Reverse Cutting

Length

	Color	cm	in	4-0	3-0	2-0	0	1	2
OS-2  19 mm 1/2 circle	Violet	14 x 14	5 x 5				RA-1001Q		
CP-2  26 mm 1/2 circle	Violet	14 x 14	5 x 5					RA-1004Q	RA-1006Q
		24 x 24	9 x 9		RA-1056Q-0			RA-1032Q	RA-1034Q
CP-1  36 mm 1/2 circle	Violet	36 x 36	14 x 14				RX-2069Q		
		45 x 45	18 x 18				RX-2068Q		
CP  40 mm 1/2 circle	Violet	36 x 36	14 x 14						RX-1066Q
		45 x 45	18 x 18						RX-2066Q



*As with traditional sutures, Quill barbed suture size is determined by its OUTER diameter, it is therefore recommended that when using Quill barbed sutures, you upsize by one USP compared with traditional suture. For example, if previously using a 3-0 traditional suture, we recommend using a 2-0 Quill barbed suture.

PDO (Polydioxanone) Bi-Directional

ABSORBABLE

Taper Point

Length

	Color	cm	in	4-0	3-0	2-0	0	1	2
RB-1  17 mm 1/2 circle	Violet	7 x 7	3 x 3	RA-2070Q					
SH-1  22 mm 1/2 circle	Violet	14 x 14	5 x 5			RA-1086Q		RA-1088Q	
SH  26 mm 1/2 circle	Violet	10 x 10	4 x 4			RA-1077Q	RA-1078Q		
		14 x 14	5 x 5				RA-1079Q		
		40 x 40	16 x 16			RA-1057Q			
 CT-1 36 mm 1/2 circle	Violet	7 x 7	3 x 3				RX-1068Q		
		14 x 14	5 x 5				RA-1000Q	RA-1003Q	RA-1005Q
		24 x 24	9 x 9				RA-1029Q	RA-1031Q	RA-1033Q
		30 x 30	12 x 12					RA-1059Q	RA-1058Q
		36 x 36	14 x 14				RA-1067Q		RX-1062Q
		45 x 45	18 x 18				RA-2067Q		
 CTX 48 mm 1/2 circle	Violet	30 x 30	12 x 12				RA-1060Q		
		36 x 36	14 x 14						RA-1065Q
		45 x 45	18 x 18						RA-2065Q
 S-18 33 mm straight	Violet	10 x 10	4 x 4				RA-1085Q		

Taper Cutting

Length

	Color	cm	in	4-0	3-0	2-0	0	1	2
 STC8 50 mm straight	Violet	7 x 7	3 x 3			RA-1013Q	RA-1023Q	RA-1024Q	

*As with traditional sutures, Quill barbed suture size is determined by its OUTER diameter, it is therefore recommended that when using Quill barbed sutures, you upsize by one USP compared with traditional suture. For example, if previously using a 3-0 traditional suture, we recommend using a 2-0 Quill barbed suture.

PDO (Polydioxanone) Uni-Directional

ABSORBABLE

Taper Point

Length

	Color	cm	in	3-0	2-0	0	1	2
RB-1  17 mm 1/2 circle	Violet	15	6	VLP-2046	VLP-2022	VLP-4109		
		20	8		VLP-2020			
SH-1  22 mm 1/2 circle	Violet	15	6		VLP-2050			
		20	8			VLP-1008		
		23	9				VLP-4149	
		30	12		VLP-2002	VLP-1003		
SKI-23  23mm SKI	Violet	23	9		VLP-2041			
SH  26 mm 1/2 circle	Violet	15	6	VLP-2045	VLP-2021	VLP-2051		
		20	8	VLP-2043	VLP-2006	VLP-2008		
		23	9					VLP-4148
		30	12	VLP-2044	VLP-2033	VLP-2009		
		45	18		VLP-2039	VLP-2038		
UR-6  26 mm 5/8 circle	Violet	15	6		VLP-2052			
		30	12				VLP-2053	

*As with traditional sutures, Quill barbed suture size is determined by its OUTER diameter, it is therefore recommended that when using Quill barbed sutures, you upsize by one USP compared with traditional suture. For example, if previously using a 3-0 traditional suture, we recommend using a 2-0 Quill barbed suture.

PDO (Polydioxanone) Uni-Directional

ABSORBABLE

Taper Point 

Length

	Color	cm	in	2-0	0	1	2
 V-30 30 mm 1/2 circle	Violet	23	9	VLP-2049			
 CT-2 26 mm 1/2 circle	Violet	15	6	VLP-4136			
		20	8	VLP-4137			
		23	9		VLP-4147		
		30	12		VLP-4135		
		70	28	VLP-4139	VLP-4140		
 CT-1 36 mm 1/2 circle	Violet	15	6				VLP-2030
		20	8		VLP-1001	VLP-2028	VLP-2016
		30	12	VLP-2042	VLP-1002	VLP-2015	VLP-2032
		45	18	VLP-2040	VLP-1004	VLP-2014	VLP-2025
		70	28		VLP-2012		VLP-2024
 CT 40 mm 1/2 circle	Violet	70	28			VLP-2013	
 CTX 48 mm 1/2 circle	Violet	45	18			VLP-2034	VLP-2037
		70	28				VLP-2010
 TP-1 65 mm 1/2 circle	Violet	45	18			VLP-2036	
		60	24				VLP-2031

*As with traditional sutures, Quill barbed suture size is determined by its OUTER diameter, it is therefore recommended that when using Quill barbed sutures, you upsize by one USP compared with traditional suture. For example, if previously using a 3-0 traditional suture, we recommend using a 2-0 Quill barbed suture.

PDO (Polydioxanone) Uni-Directional

ABSORBABLE

Precision Reverse Cutting

Length

	Color	cm	in	3-0	2-0	0	1	2
P-3  13 mm 3/8 circle	Violet	20	8	VLP-3021				
		30	12	VLP-3022				
PS-2  19 mm 3/8 circle	Violet	20	8	VLP-3026	VLP-3027			
		30	12	VLP-4132	VLP-3028			
		45	18	VLP-2029	VLP-2001			
PS-1  24 mm 3/8 circle	Violet	30	12		VLP-3029			
		45	18		VLP-2026			
		60	24		VLP-2027			
PS  26 mm 3/8 circle	Violet	20	8			VLP-1007		
		30	12			VLP-2005		

Reverse Cutting

Length

	Color	cm	in	3-0	2-0	0	1	2
CP-1  36 mm 1/2 circle	Violet	20	8			VLP-1005		
		30	12			VLP-1006		
		45	18			VLP-4138	VLP-4144	
		60	24			VLP-4143		
		70	28				VLP-4146	
CP  40 mm 1/2 circle	Violet	45	18					VLP-2019
		70	28					VLP-2011
OS-6  36 mm 1/2 circle	Violet	70	28			VLP-4142		
OS-8  40 mm 1/2 circle	Violet	60	24				VLP-4134	

*As with traditional sutures, Quill barbed suture size is determined by its OUTER diameter, it is therefore recommended that when using Quill barbed sutures, you upsize by one USP compared with traditional suture. For example, if previously using a 3-0 traditional suture, we recommend using a 2-0 Quill barbed suture.